

# Spatial gradient contour maps for **Artist** systems


This contour map shows the spatial gradient for **Artist** system (which is a 1.5T system with 70cm patient bore and **DVw** magnet.)

The map covers a range of  $\pm 2$  meters from isocenter along the magnet's axis, and  $\pm 1.6$  meters across the magnet.


Minimum entry spatial gradient:

To reach the center of the magnet requires passing through at least a SG=**291** G/cm.

The contour locations and patient bore dimensions are accurate to isocenter. However, the enclosure's exterior shape is approximate, due to differences between enclosure designs.



# Detail contour maps for Artist systems



**D** can be read from the map, the table, or the plot. At 300 G/cm, **D** is 26.8 cm.

Spatial Gradient, G/cm	Distance <b>D</b> , cm
<291	35*
300	26.8
350	15.9
400	10.9
450	7.5
500	5.3
550	3.5
600	2.1
650	1.0
700	0.1
≥710	0

\* This is the minimum entry SG, where **D** reaches half of the bore diameter.

This contour map zooms in on the area near the mouth of the magnet. The map shows how far a given contour reaches pasting the bore wall. This is distance **D** – it can be read from the contour map, from the table or from the plot. The enclosure shape is accurate in this plot.

An example – marked in red in the plot, the table, and the graph – shows that for SG of 300 G/cm, distance **D** is 26.8 cm. This leaves an “opening” of 16.4 cm.

